

JumpStart® Kangaroo®

AFO trimline, PF block, DF free



CASCADE
daFO®

The JumpStart Kangaroo was developed to help young patients with fully-correctable pronation or supination and excess plantarflexion re-align their base of support to improve stability in standing and walking. Use of the Kangaroo encourages more practice in dorsiflexed positions for patients who present with excess toe walking or knee hyperextension. Full wrap-around flexible support encompasses the foot and ankle for pronation and supination control, and the well-defined plantar surface provides proprioceptive input. Proximal, posterior AFO trimline blocks excess plantarflexion, encouraging knee flexion and tibial movement over the foot. Solid outer foundation adds extra arch support, heel stability and strength for block to excess plantarflexion.



Key Features

Open anterior

allows access to free dorsiflexion.

Solid outer foundation and posterior leg portion

add extra arch support, heel stability, and strength to block excess plantarflexion.

Full wrap inner liner

made from soft, thin polyethylene plastic allows full alignment control of the heel, mid-foot, and forefoot.

Choice of two outer foundation materials

Polyethylene for moderate flexibility or co-poly for firm support.

Optional Super Heel Seating Strap

provides additional plantarflexion block and heel control.

Prefabricated

sized by measurement using the Fast Fit sizing jig.

Widths: Wide and Narrow

Sizes: 4.00 – 9 inches (0.25 inch increments)

Designed for patients who:

- Exhibit a fully correctable foot position and fit within the product size range.
- Present with moderate pronation or supination and associated gait instability.
- Require block to excess plantarflexion presented as strong toe walking or knee hyperextension.
- May benefit from the sensory feedback provided by a fully contoured foot bed.
- Will benefit from full access to dorsiflexion.



No casting needed!