

The DAFO 9 Softy is designed for limited daytime use and/or as part of a night-stretching program. The stretching straps of the articulated proximal section allow positional adjustment in order to help gain or maintain ankle range. The DAFO 9 Softy provides excellent hindfoot and forefoot control, promoting proper foot positioning. The soft foam liner is ideal for patients with challenging foot anatomy, difficult-to-correct feet, or those prone to pressure points and irritation.



Key Features

Custom-fabricated from a cast or scan

Magnetic fastener

for easy donning and doffing

Extended toe shelf trimline

allows even pull of entire foot-bed

Designed for patients who:

- lack range of motion at the ankle due to tight calf muscles and Achilles tendon.
- would benefit from passive, adjustable stretching when resting or sleeping to improve range of motion at the ankle.



Fastener attached



Fastener unattached